



HEINE HSL 150® Hand-held Slit Lamp

For the examination of the anterior segment



Compact and light-weight, the portable HSL 150 is ideal for situations where a full-size slit lamp is not available or convenient. For the examination of the anterior segment in Human- and Veterinary medicine.

- ∴ **Slit size adjustable** from 0.2 mm x 10 mm to 4 mm x 14 mm.
- ∴ **6x magnification.**
- ∴ **Blue interference filter (FITC)**, can be selected for corneal examination.
- ∴ **Multi-coated optics** for maximum light transmission.
- ∴ **2.5V or 3.5V XHL Xenon Halogen Technology** for bright, white light comparable with the brightness of a classic slit lamp.
- ∴ **The HSL 150 instrument head weighs just 70 g**, the lightest slit lamp of its kind.
- ∴ **Available with HEINE's exclusive AV connector or TL "twist lock" style.** AV features a positive lock, TL offers compatibility with existing equipment.

HSL 150 for BETA or BETA SLIM handles	AV	TL
HSL 150, head only	C-02.14.602	C-02.14.602 TL
XHL Xenon Halogen spare bulb	X-02.88.099	X-02.88.099
Replacement filter attachment for HSL 100 / 150		C-00.14.605

HEINE alpha+ HSL 150® Hand-held Slit Lamp

Superior full-size performance in a convenient size



For the examination of the anterior segment in Human- and Veterinary medicine. Basic technical specification similar to HSL 150 (see above), but can only be used with the alpha+ handle system.

alpha+ HSL 150 for alpha+ handles	2.5 V	3.5 V
alpha+ HSL 150, without handle	D-01.81.511	D-02.81.511
alpha+ HSL 150 with alpha+ handle		
alpha+ HSL 150 with alpha+ battery handle	D-302.10.051	
alpha+ HSL 150 with alpha+ rechargeable handle and NT 200 B charger		D-302.20.444
XHL Xenon Halogen spare bulb	X-01.88.098	X-02.88.099
Replacement filter attachment for alpha+ HSL 150		D-00.81.502

HEINE HSL 10x® Loupe Attachment

Loupe attachment with 10x magnification



For HSL 100, HSL 150 and alpha+ HSL 150. With precision magnification optics for crisp, clear images.

- ∴ **10x magnification** when combined with HSL 100, HSL 150 or alpha+ HSL 150.
- ∴ **Increased working distance** from examiner to patient.
- ∴ **Large field of view.**
- ∴ **With an additional detachable TP- elastomer Eyepiece** to steady the instrument on the brow and avoid stray light. Can be replaced with a soft TP- elastomer cap, suitable for spectacle-wearers.

HSL 10x Loupe Attachment	C-00.14.606
--------------------------	-------------