

InFocus

Look closer.
See further.

BX 900 LED

for professional photography

InFocus August 2014

Special topic about features and devices or services from Haag-Streit Diagnostics.

Published by
Haag-Streit AG
Marketing & Sales
Gartenstadtstrasse 10
3098 Koeniz
Switzerland

New LED powered BX 900 Slit lamp

BX 900 LED



It is our great pleasure to announce the introduction of the BX 900 LED. After the clinical slit lamps BQ 900, BP 900 and BM 900 now also the photographer slit lamp BX 900 benefits from the advantages of the LED technology.

BX 900 LED

- Ergonomic one-hand control for slit and background illumination
- Most homogeneous slit for best examination
- Illumination lasts a lifetime – no more change of light bulbs
- Energy saving
- EyeSuite Imaging

Advanced usability – improved illumination

BX 900 LED



Ergonomic one-hand control

The illumination head delivers independently light for slit and background which can be continuously variable controlled by a single control panel. This allows adjusting slit and background illumination with one hand which benefits the task of the photographer.

Standard LED illumination

The LED slit illumination is meanwhile well established and delivers the most homogeneous slit ever, thus providing more details for diagnosis. The illumination for the photographs taken with the BX 900 remains unchanged – it is still delivered from our very powerful flash unit.

EyeSuite Imaging


From the latest release on the automated import of images in various formats is supported* and can be activated with a dongle license code. This allows that now also the BX 900 can be integrated in the standard Haag-Streit EyeSuite software.



Ergonomic one-hand control



Standard LED illumination



EyeSuite Imaging

* Images in formats BMP, TIFF, PNG, JPEG in size <18 MPXventional systems