



# HandyRef-K / HandyRef

HANDHELD REF/KERATOMETER / HANDHELD REFRACTOMETER

**Anytime-anywhere** ✓ handheld measurement

Improved usability by **valiant contraptions** ✓

**Accurate measurement** ✓ with quality functions

Interface **enhancement** ✓

✓ **Superior** optional accessories for **even easier** ✓ use



Children chart is optional.

## HandyRef-K / HandyRef Specifications

<b>Main body</b>	
Display	3.5-inch color LCD
Interface	USB: 1 port, Wireless LAN (WLAN): 1ch (only for regions approved by Radio Law)
<b>Power Specification</b>	
Battery pack	Lithium-ion battery (7.4 V 1850 mAh)
Station feed	DC 9 V 2 A (maximum)
<b>Dimensions</b>	
	205 (W) x 181 (D) x 224 (H) mm (including occluders) 8.1 (W) x 7.1 (D) x 8.8 (H)" (including occluders)
<b>Station</b>	
Printer	Thermal line printer with auto cutter (printer-equipped model only)
Interface	USB: 1 port, LAN: 1 port, RS-232C: 1 port (printer-equipped model only)
Power supply	AC 100 to 240 V, 50 / 60 Hz
Power consumption	60 VA
<b>Dimensions / Mass</b>	
	224 (W) x 283 (D) x 147 (H) mm / 3 kg (printer-equipped model), 2.5 kg (model without printer) 8.8 (W) x 11.1 (D) x 5.8 (H)" / 6.6 lbs. (printer-equipped model), 5.5 lbs. (model without printer)
<b>Standard accessories</b>	
	Occluder (2 units), Neck strap, Printer paper (3 rolls / printer-equipped model only), Power cord, Connection cable, Lithium-ion battery, Dust cover, Spherical model eye (contact lens holder integral type)
<b>Optional accessories</b>	
	Children chart, Carrying case, Carrying case with stand, Barcode scanner, Magnetic card reader, EyeCa-RW2, Eye Care card, Communication cable, Lithium-ion battery, USB flash drive

CE Mark for this product is pending.

This product is not cleared by the FDA for distribution in the United States.

Specifications may vary depending on circumstances in each country.

Specifications and design are subject to change without notice.



**HEAD OFFICE**  
34-14 Maehama, Hiroishi  
Gamagori, Aichi 443-0038, Japan  
Telephone: +81-533-67-6611  
Facsimile : +81-533-67-6610  
URL : <http://www.nidek.co.jp>  
[Manufacturer]

**TOKYO OFFICE  
(International Div.)**  
3F Sumitomo Fudosan Hongo Bldg.,  
3-22-5 Hongo, Bunkyo-ku, Tokyo  
113-0033, Japan  
Telephone: +81-3-5844-2641  
Facsimile : +81-3-5844-2642  
URL : <http://www.nidek.com>

**NIDEK INC.**  
47651 Westinghouse Drive  
Fremont, CA 94539, U.S.A.  
Telephone: +1-510-226-5700  
: +1-800-223-9044 (US only)  
Facsimile : +1-510-226-5750  
URL : <http://usa.nidek.com>

**NIDEK S.A.**  
Europarc  
13, rue Auguste Perret  
94042 Créteil, France  
Telephone: +33-1-49 80 97 97  
Facsimile : +33-1-49 80 32 08  
URL : <http://www.nidek.fr>

**NIDEK TECHNOLOGIES Srl**  
Via dell'Artigianato, 6 / A  
35020 Albignasego (Padova), Italy  
Telephone: +39 049 8629200 / 8626399  
Facsimile : +39 049 8626824  
URL : <http://www.nidektechnologies.it>

